## Self-Cured Dental Universal Adhesive REF AD03 TOKUYAMA UNIVERSAL BOND

### ENGLISH

Read all information, precautions and notes before using.

#### ■PRODUCT DESCRIPTION AND GENERAL INFORMATION

- TOKUYAMA UNIVERSAL BOND is a two-component self-cured dental adhesive system for both direct and indirect restorations that can be used with Self-Etch, Selective-etch, etch and Total-etch techniques. As a universal adhesive, TOKUYAMA UNIVERSAL BOND has been designed to be fully compatible with light-cured, self-cured and dual-cured composite materials.
- TOKUYAMA UNIVERSAL BOND improves the bond strength of polymerizable resin material (adhesive resin cement, acrylic resin and composite resin) to indirect restorative materials such as glass-ceramics (porcelain), oxide-ceramics (zirconia and alumina), metals (precious and non-precious) and resin materials including inorganic filler.
- TOKUYAMA UNIVERSAL BOND contains Phosphoric acid monomer, Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylen glycol dimethacrylate (TEGDMA), 2-Hydroxyethyl methacrylate (HEMA), MTU-6 (Thiouacil monomer), Silane coupling agent, Peroxide, Borate catalyst, Acetone, Isopropanol and purified water.

#### ■INDICATIONS

- Direct anterior and posterior restorations with light-curing, dual-curing, and self-curing composite materials
- Intraoral repair of composite restorations, porcelain fused to metal, metal, and all-ceramic restorations without an additional primer
- Cementation of indirect restorations and veneers when combined with light-cure, dual-cure, and self-curing resin cement
- Bonding of core build-ups made of core build-up materials
- Bonding of denture resin to metal base, clasp or attachment
- Repair of denture with metal base, clasp or attachment
- Bonding of opaque resin to a metal base in the fabrication of resin-faced stainless steel crowns

#### ■CONTRAINDICATIONS

TOKUYAMA UNIVERSAL BOND contains methacrylic monomers, organic solvents and acids. DO NOT use TOKUYAMA UNIVERSAL BOND for patients allergic to or hypersensitive to methacrylic monomers, related monomers, organic solvents and acids.

#### ■PRECAUTIONS

- DO NOT use TOKUYAMA UNIVERSAL BOND for any purpose other than those listed in these instructions. Use TOKUYAMA UNIVERSAL BOND only as directed herein.
- TOKUYAMA UNIVERSAL BOND is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- DO NOT use TOKUYAMA UNIVERSAL BOND if the safety seals are broken or appear to have been tampered with.
- If TOKUYAMA UNIVERSAL BOND causes an allergic reaction or oversensitivity, discontinue its use immediately.
- Use examination gloves (plastic, vinyl or latex) at all times when handling TOKUYAMA UNIVERSAL BOND to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If TOKUYAMA UNIVERSAL BOND comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- Avoid contact of TOKUYAMA UNIVERSAL BOND with eyes, mucosal membrane, skin and clothing.
- If TOKUYAMA UNIVERSAL BOND comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
- If TOKUYAMA UNIVERSAL BOND comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the completing the restoration. Affected areas may whiten from protein coagulation; however, whitening should disappear within 24 hours. If whitening does not disappear within 24 hours, immediately contact a physician, and the patient should be so advised.
- If TOKUYAMA UNIVERSAL BOND comes into contact with the skin or clothing, immediately saturate the area with an alcohol-soaked cotton swab or gauze and flush with water.

- Instruct the patient to rinse his mouth immediately after treatment.
- TOKUYAMA UNIVERSAL BOND should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- To avoid the unintentional ingestion of TOKUYAMA UNIVERSAL BOND, do not leave it unattended within the reach of patients and children.
- DO NOT expose TOKUYAMA UNIVERSAL BOND or its vapor to open flame because TOKUYAMA UNIVERSAL BOND is flammable.
- To avoid cross infection or lowering bond strength, DO NOT reuse the disposable applicator. Clean the mixing well thoroughly with alcohol after each use.
- DO NOT use TOKUYAMA UNIVERSAL BOND to inject into the cavity being in close proximity to the pulp. Pulp protection with glass ionomer lining or calcium hydroxide is recommended. DO NOT use EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of TOKUYAMA UNIVERSAL BOND.
- The substances listed below inhibit adhesion and must be removed from the adhesive surface by mechanical cleaning or alcohol soaking. Clean the restoration after adjustment if necessary.
- DO NOT use mixing wells supplied by other manufacturers.
- Use the mixing well or disposable mixing well provided when dispensing and mixing of TOKUYAMA UNIVERSAL BOND. If TOKUYAMA UNIVERSAL BOND is not dispensed into the rubber tip it may be difficult to remove the cured TOKUYAMA UNIVERSAL BOND from mixing well. The cured TOKUYAMA UNIVERSAL BOND can be removed by inverting the rubber.
- DO NOT mix TOKUYAMA UNIVERSAL BOND with other brands of primers or sealers.

- To avoid splashing the liquid, the bottle of TOKUYAMA UNIVERSAL BOND should not be opened immediately after taking out of the refrigerator. TOKUYAMA UNIVERSAL BOND should be brought to room temperature before using. After removing from the refrigerator, allow the product to sit for 20 minutes or until it reaches room temperature.
- PRECAUTIONS FOR MEDICATIONS AND MATERIALS**
  - Some materials and medications (hemostatic material) inhibit setting/adhesion of TOKUYAMA UNIVERSAL BOND.
    - DO NOT USE products that contain:
      - eugenol,
      - aluminum chloride,
      - ferriic sulfate,
      - calcium sulfate,
      - aluminum sulfate,
      - fluamine silver fluoride [molecular formula: Ag(NH<sub>4</sub>)F<sub>2</sub>]
  - TOKUYAMA UNIVERSAL BOND does not adhere to tooth structure immediately after applications of materials that contain:
    - hydrogen peroxide (toxold),
    - sodium hypochlorite
- We advise avoiding the use of these materials on the day of cementation, however these materials can be used for root canal treatment as the inhibition caused by these materials usually evaporates within five days.

#### ■STORAGE

- Store TOKUYAMA UNIVERSAL BOND under refrigeration at temperatures between 0 to 10°C (32 to 50°F).
- Keep away from heat, direct sunlight, sparks or open flames.
- DO NOT use after the indicated expiration date on the bottle/package.

#### ■DISPOSAL

Unused TOKUYAMA UNIVERSAL BOND is to be absorbed into an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

#### ■CLINICAL APPLICATION FOR DIRECT RESTORATIONS AND INTRAORAL REPAIR OF RESTORATIONS WITH COMPOSITE RESIN

#### 1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

#### 2. Isolation

The rubber dam is the preferred method of isolation.

#### 3. Cavity preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (classes I, II) because bevels and chamfers assist in easing demarcations between the cavity margins and the restoration, thereby enhancing both esthetics and retention.

- If the adherent surface involves uncut enamel, apply etching agent to the uncut enamel. See section "Selective-enamel-etch technique".

#### 4. Drying

Dry the cavity by using a blotting technique or an air syringe technique.

- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
- The substances listed below, which inhibit curing of TOKUYAMA UNIVERSAL BOND , must be removed from the tooth surface by thoroughly cleaning the tooth surface with alcohol, citric acid, or the application of phosphoric acid for 2 to 3 seconds before application:

- Oil mist from handpiece
- Plaque, tartar, saliva, blood and exudate fluids
- If necessary, protect the surface not being bonded with metal separator, paraffin film, Vaseline, etc.

#### 5. Etching of the tooth structure (optional)

Clinically adequate adhesive bonds are achieved by the application of the adhesive. Selectively etching the enamel or using the "total etching procedure" can achieve higher adhesive values.

- Selective-enamel-etch technique
  - Apply phosphoric acid gel only to the enamel surrounding the margin of prepared cavity and leave phosphoric acid gel in place for 10-15 seconds. Rinse the etched surface thoroughly (at least 15 seconds) with water, and then dry with mild air.
- Total etch technique
  - Apply phosphoric acid gel onto the prepared enamel first, and then the dentin. The etchant should be left to react on the enamel and dentin for 10-15 seconds. Rinse thoroughly (at least 15 seconds) with water and dry with mild air or with cotton pellets, do not overdry.

#### 6. Preparing the adhesive surface during repair

- Ceramics other than Porcelain, and Composite materials
  - Bevel the margin if necessary. Roughen the area using air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.
- Metal restorations
  - Roughen the area using air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.
- Porcelain
  - Bevel the margin if necessary. DO NOT touch the surface; Polish the area to prepare a clean and fresh surface for adhesion. Thoroughly rinse and dry the surface.
  - When using air abrasion, use at least 30 to 50 micron grid of alumina powder.
  - The substances listed below inhibit adhesion and must be removed from the adhesive surface by mechanical cleaning or alcohol soaking.
  - Oil mist from handpiece
  - Plaque, tartar, saliva, blood and exudate fluids
  - If necessary, protect the surface not being bonded with metal separator, paraffin film, Vaseline, etc.

diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

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#### 2. Dispensing of TOKUYAMA UNIVERSAL BOND

Allow TOKUYAMA UNIVERSAL BOND to reach room temperature prior to use. Open the bottle by removing the cap and dispense one drop each of TOKUYAMA UNIVERSAL BOND A and B into the mixing well or disposable mixing well. Mix thoroughly with a disposable applicator.

- Complete the application within 1 minute of dispensing when using the mixing well since TOKUYAMA UNIVERSAL BOND contains volatile solvents. After mixing, the color of TOKUYAMA UNIVERSAL BOND changes gradually based on the reaction of the catalyst of TOKUYAMA UNIVERSAL BOND, this is normal and does not indicate a problem with TOKUYAMA UNIVERSAL BOND.
- Complete the application within 3 minutes of dispensing when using the disposable mixing well.

- Be sure to hold the bottles of TOKUYAMA UNIVERSAL BOND A and B vertically when dispensing. Do not hold the bottles diagonally or horizontally while dispensing; this could cause the TOKUYAMA UNIVERSAL BOND flow back into the nozzle and become contaminated.

- Wipe off excess TOKUYAMA UNIVERSAL BOND from the tip of the nozzle prior to closing.
- Immediately after dispensing, tightly replace the bottle caps of TOKUYAMA UNIVERSAL BOND A and B.

#### 3. Application of TOKUYAMA UNIVERSAL BOND

Using the disposable applicator, apply the mixed TOKUYAMA UNIVERSAL BOND on the surface to be bonded. Ensure TOKUYAMA UNIVERSAL BOND covers all surfaces where the cement will be applied.

- DO NOT over apply. One application of TOKUYAMA UNIVERSAL BOND is sufficient. Too thick of a layer may compromise the fit of the restoration.
- DO NOT allow splashing of TOKUYAMA UNIVERSAL BOND to come in contact with the surface to be bonded. If the treated surface is contaminated, thoroughly rinse with water, dry and re-apply fresh TOKUYAMA UNIVERSAL BOND.
- DO NOT wait too long after the application of TOKUYAMA UNIVERSAL BOND. Air dry within 30 seconds after application, to ensure the appropriate film thickness.

#### 4. Air Dry

Using an oil-free air/water syringe, apply mild air to the surface. Use a vacuum aspirator to prevent spatter of the TOKUYAMA UNIVERSAL BOND.

- DO NOT use the dry surface after TOKUYAMA UNIVERSAL BOND is applied.

- Cementing**

Cement with a composite resin cement according to its manufacturer's instructions.

#### ■CLINICAL APPLICATION FOR BONDING OF ACRYLIC RESIN TO METAL BASE, CLASP OR ATTACHMENT

1. **Preparing the adhesive surface**

- Roughen the area by sandblasting (0.1 to 0.2MPa), or by grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

- Dispensing, Application and Air Dry**
  - Use the same procedures as described in CLINICAL APPLICATION FOR CEMENTATION OF INDIRECT RESTORATIONS.

- Application, curing, trimming and finishing of acrylic resin**
  - Proceed with the procedure according to the manufacturer's instructions for the use of polymerizable acrylic resin.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use. Specifications are subject to change without notice. When the product specifications change, the instructions and precautions may also change.

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