AKZENTA

Disposable Nitrile Examination Gloves

GB How to use the gloves
SV Instruktioner för användning
DK Brugsanvisning
NO Instruksjon for bruk
FI Käyttöohjeet
DE Gebrauchsanleitung
PL Instrukcja użytkowania
FR Mode d'emploi
IT Istruzioni d'uso
ES Instrucciones de uso
NL Gebruiksaanwijzing
BG Инструкции за употреба

CS Návod k použití
HR Upute za uporabu
ET Kasutusjuhend
EL Οδηγιο χρησης
LV Lietošanas instrukcija
LT Naudojimo instrukcija
PT Instruções de utilização
RO Instrucțiuni de utilizare
SK Návod na použitie
SL Navodila za uporabo
HU Használati utasítás

Rev. 20092024

 Cat III Reg (EU) 2016/425
Class I MDR (EU) 2017/745

For single use only

EU DoC can be downloaded at : www.akzenta.com

Manufactured By
Akzenta International SA Via Giuseppe Motta, 24, CH-6830 Chiasso. Phone: +41 (0) 91 921 14 92 www.akzenta.com
Importer
IQX GmbH Kudlichstraße 37, 4020 Linz, Austria.

Art. No.	Model	Size	Standard Package
GLOBY53-AKZ-TL	Nitrile Supreme Style Edition	XS/S/M/L	1/10

AQL 1,5 EN 455 - 1, 2, 3, 4

These gloves are tester and approved by Notified Body: 2777 module B
SATRA Technology Europe Ltd. Bracetown Business Park, Clonee, Dublin 15, Dublin, Ireland.

Notified Body responsible for Module D Conformity Assessment:
SATRA Technology Europe Ltd. Notified Body No. 2777 Bracetown Business Park, Clonee, Dublin 15, Dublin, Ireland.

EN ISO 374-1:2016+A1:2018			
Measured breakthrough time (min)		Permeation performance level	
>10		1	
>30		2	
>60		3	
>120		4	
>240		5	
>480		6	

EN ISO 374-1:2016+A1:2018 Type C		EN ISO 374-4:2019	EN ISO 374-5:2016
Permeation levels are based on breakthrough times		Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical	The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen
Breakthrough Performance Level		Resistance to chemical degradation (%)	
GLOBY53-AKZ-TL		-25.2	Protection against bacteria and fungi - Pass
40% Sodium Hydroxide (K)		6.5	Protection against viruses - Pass
30% Hydrogen Peroxide (P)		17.9	
37% Formaldehyde (T)			

GB

EN ISO 374-1:2016+A1:2018	EN ISO 374-1:2016+A1:2018	EN ISO 374-1:2016+A1:2018
Type A	Type B	Type C
		
UVWXYZ	XYZ	X
Type A - At least Breakthrough Performance Level Class 2 (more than 30 minutes) against at least 6 chemicals on the list. * Fig. 1	Type B - At least Breakthrough Performance Level Class 2 (more than 30 minutes) against at least 3 chemicals on the list. * Fig. 1	Type C - At least Breakthrough Performance Level Class 1 (more than 10 minutes) against at least 1 chemicals on the list. * Fig. 1

The EU Declaration of Conformity is accessible at www.akzenta.com

Product description
Non-Sterile Powder Free Nitrile Examination Gloves are made from 100% synthetic rubber, ambidextrous, coloured and for single use. Gloves have high strength and flexibility, while at the same time they are soft and well retain the original shape. Gloves are ambidextrous with a comfortable rim that ensures easy donning, ideal fit, secure fixing (the glove does not slip) and the rapid removal of gloves at the end of the manipulation.
These gloves also contribute to moisture the skin of the user and have been tested according to EN ISO 21420:2020 standard.
These gloves are to protect user's hands against microorganism (bacteria, fungus and viruses) and certain chemical risks. Testing and marking on the gloves are in accordance with Regulation (EU) 2016/425 as well as applicable harmonized European Standards. The gloves are to be used strictly for the intended applications. Potential consumers of these gloves may be personnel from medical and preventive institutions, private clinics, dentist, nurses, cosmetologist, lab personnel, food industry personnel and forensic personnel.
Additional information on product constituents will be provided upon request.

This product is landfill biodegradable. The actual biodegradation rates

will vary depending on the landfill conditions and the biological activity of microorganism surrounding the nitrile gloves.
These gloves have been tested per ASTM 5526 testing and validates up to 75% biodegradation within 776 days.

Storage conditions
Do not store Non Sterile Powder Free Nitrile Examination Gloves where temperature may rise above 104°F (40°C). Store in cool, dry and well ventilated area. Opened boxes of Non Sterile Powder Free Nitrile Examination Gloves should be shielded from exposure to direct sunlight or fluorescent lighting to prevent discoloration. Improper storage of Non Sterile Powder Free Nitrile Examination Gloves will result in decreased shelf life and compromised efficiency.

Use
These gloves are designed to protect user's hands against certain chemical risks. Testing and marking on the gloves are in accordance with Regulation (EU) 2016/425 as well as applicable harmonized European Standards. Please ensure the gloves are used strictly for the intended applications.

Product instructions for use

- Before usage, inspect the gloves for any defect or imperfections. Usage is for single use only. If re-used, the risk of contamination and infection increases due to improper cleaning processes and increased risk of holes and tear during re-use due to weakening of gloves by cleaning processes.
- Sizing: Select the right size of glove for your hands.
- Donning: Hold the glove by the bead with one hand. Align the glove thumb with your other hand thumb and slide your hand into the glove, one finger into each glove finger. Pull by the glove palm to get a good fit. Don the other glove by the same procedure.
- Inspection: punctures or tears may occur after donning. Inspect the glove after donning and stop to use them immediately if a damage is found.
- Doffing: Hold the glove bead and pull toward the finger until the glove comes off.
- Disposal: For the proper disposal of all used gloves, please follow your institution's policies for disposal.

Caution

Primary material – Nitrile rubber. Does not contain natural rubber latex.

Components used in glove manufacturing may cause allergic reactions in some users.

Warning

EN ISO 374-4:2019 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm- where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture.

It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

Before usage, inspect the gloves for any defect or imperfections. Discard when the gloves are damaged, FOR SINGLE USE ONLY.

The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

Cleaning

Not applicable as the glove is single use in intended manner.

Shelf Life

The shelf life of the product is 5 years from the date of manufacture.

Reporting

In the event any serious incident occurs in relation to the use of this device, please report the incident to the manufacturer or its authorized representative and the competent authority of the Member State in which the user is established.

Dispose of the gloves and all biologically contaminated matter in the appropriate container. Product biodegradation applicable in landfill disposal only.

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