


Endo Instrument Holders



Conformity / Regulatory Information

Generic Device Group	Endo Instrument Holders	
Model	1180	
Regulation	Medical Device Regulation (EU) 2017/745	
Classification and Rule	Class I, Rule 1	
Basic UDI-DI	++D125ENDOINSTHOLDERSGF	
Nomenclature	EMDN	Q010199
	GMDN	64813
	UMDNS	16354
Intended Use	A device intended to be used during an endodontic procedure to facilitate the access, handling, temporary storage and convenient selection of endodontic instruments. The holder is intended to be used with disposable spongy foam inserts for holding and cleaning delicate and sharp rotary endpieces.	
Indications of use	- For professional dental application only	

Product Specification

Dimensions	50 x 50 x 25 mm (L x W x H)	
Features & Accessories	<ul style="list-style-type: none"> - Built in endo file ruler (up to 34 mm) - For use with: 801 Square Endo Cushions (50 pcs/pack) 	
Packaging Description	Primary	1 pc/pack
	Secondary	100 packs/case
Packaging Languages	EN, PL, DE, FR, IT, ES, PT, FI, NO, SV, DA, RO, CS, EL, HU, LV, LT, SK, JA, ZH	
Composition	Stainless Steel	
Traceability	UDI / LOT	
Storage Conditions	Do not store in direct sunlight, store in a cool dry place	
Reprocessing	Autoclavable to 134°C	

Manufacturer Information

Manufacturer	Premium Plus (Dongguan) Ltd.	Manufacturer SRN	CN-MF-000038259
EU Representative	Premium Plus Poland sp. z o.o.	EU Representative SRN	PL-AR-000030110
UK Representative	Premium Plus UK Ltd.	Country of Manufacture	China

Premium Plus (Dongguan) Ltd.
No.1 Industrial Area, Tutang,
Changping, Dongguan,
Guangdong 523581
China
Tel: 86-769-83397277
Email: info@premiumpluschina.com

Premium Plus Poland sp. z o.o.
ul. Bukowska 27
62-081 Wysogotowo
Poland
Tel: 48-61-880-1094
Email: info@premiumpluspl.com
Web: www.premiumpluspoland.com

Premium Plus UK Ltd.
2, Knighton Heath Industrial Estate
847 Ringwood Road, Bournemouth,
Dorset BH11 8NE
U.K.
Tel: 44-1202-610101
Email: info@premiumplusuk.com
Web: www.premiumplusuk.com

Premium Plus Medical Supplies Ltd.
1001, Yuen Long Trading Centre
No.33, Wang Yip Street West
Yuen Long, N.T.
Hong Kong
Tel: 852-2443-1219
Email: hk@premiumplushk.com
Web: www.premiumplushk.com

Premium Plus Dental Supplies Inc.
155 Liberty Avenue
Brooklyn NY 11212
U.S.A.
Tel: 1-718-366-8686
Email: sales@premiumplususa.com
Web: www.premiumplususa.com

Premium Plus Japan Co., Ltd.
5th Floor, 1-1-31, Nishimidorigaoka,
Toyonakashi, Osaka 560-0005
Japan
Tel: 81-6-6845-0066
Email:
premiumplus.info@premiumplus.jp
Web: www.premiumplus.jp