

ProphyCare® HAp

ProphyCare HAp is a fluoride and paraben free, one-step prophy paste that effectively removes stains while helping to reduce hypersensitivity.



**Fluoride and
Paraben free**



Smooth & Creamy

The appealing mild mint taste and creamy texture is pleasant for patients of all ages.



One Step Procedure

ProphyCare HAp is a one-step prophy paste offering effective stain removal in the beginning of the treatment and giving the paste its final gentle polishing properties in the end.




Remineralizing

Studies* have shown that products containing hydroxyapatite have better remineralization properties than those of traditional fluoride products.


*Enamel and dentine remineralization by nano-hydroxyapatite toothpastes, JOD, Volume 39, Issue 6, June 2011, Pages 430-437, study by Peter Tshoppe, Daniela L.Zandim, Peter Martus, Andrej M.Kielbassa

Desensitisation

ProphyCare HAp has a composition containing hydroxyapatite that effectively blocks the open dentin tubules on an exposed dentin surface, helping to reduce hypersensitivity. (5KX magnification. Data on file)



Untreated dentin surface



Dentin surface after polishing with ProphyCare HAp

**“Contains hydroxyapatite
that contributes to relief from
hypersensitivity”**



For more information scan the QR-code with your mobile phone or visit www.directadental.com

ProphyCare® Blue 250
 Coarse



ProphyCare® Blue 250
 Art: 731113 tube, 60ml
 Art: 690123 single dose, 144 x 2g

**“The perfect starter,
 effective removal
 of stubborn stains”**

ProphyCare® Green 170
 Medium



ProphyCare® Green 170
 Art: 731114 tube, 60ml
 Art: 690122 single dose, 144 x 2g

**“Less abrasive than Blue,
 but still effective”**

ProphyCare® Red 120
 Fine



ProphyCare® Red 120
 Art: 731112 tube, 60ml
 Art: 690121 single dose, 144 x 2g

**“Gentle cleaning, tissue
 friendly stain remover”**

ProphyCare® Yellow 40
 Extra Fine



ProphyCare® Yellow 40
 Art: 731111 tube, 60ml
 Art: 690120 single dose, 144 x 2g

**“Ideal for
 final polishing”**

ProphyCare® Paraben free



**ProphyCare® Blue 250
 Paraben free**
 Art: 731116 tube, 60ml
 Art: 690103 single dose, 144 x 2g



**ProphyCare® Green 170
 Paraben free**
 Art: 731119 tube, 60ml
 Art: 690102 single dose, 144 x 2g



**ProphyCare® Red 120
 Paraben free**
 Art: 731115 tube, 60ml
 Art: 690101 single dose, 144 x 2g



**ProphyCare® Yellow 40
 Paraben free**
 Art: 731118 tube, 60ml
 Art: 690100 single dose, 144 x 2g

ProphyCare® PRO

ProphyCare® PRO
 One-step prophy paste

**“Integrated abrasion
 variability, combining plaque
 removal and polishing”**

**ProphyCare® PRO
 Lilac Universal**
 Art: 741111 tube, 60ml
 Art: 690125 single dose, 144 x 2g



ProphyCare® HAp

ProphyCare® HAp
 One-step prophy paste

**“Contains hydroxyapatite
 that contributes to relief from
 hypersensitivity”**

ProphyCare® HAp
 Art: 731117 tube, 60g
 Art: 690108 single dose, 144 x 2g

